




PT100 Thin Probe

Probe designed to measure temperature
in incubator inserts with lids (IVF)



Non-contractual photo

Réf. 13357

Introduction

The PT100 Fine probe measures the temperature in incubator inserts with lids (IVF) up to 600°C. It is particularly suited to the needs of industrial and pharmaceutical applications.

Easy to use, the probe connects directly to the Nano SPY U recorder (provide the connection pack ref. 12617) or to the LoRa SPY U. It does not require the use of a converter or any external power supply.

Compatible with version 4.0 of the Nano SPY U.


Technical features

Measurement range	From -50°C to +600°C of the sensitive element
Sensor class	A - CEI 60751
Sensor dimensions	Ø1.5x120mm
Cable type and dimensions	3-wire pvc/pvc cable - 3m
Cable operating temperature	From -50°C to +600°C

Connection

Case 1: Connection diagram with a LoRa® SPY U


- White wire of the PT100 IVF sensor: Connect the extremity to connector 5 of the LoRa® SPY U.
- 1st Red wire of the PT100 IVF probe: Connect the extremity to the connector 3 of the LoRa® SPY U.
- 2nd Red wire of the PT100 IVF probe: Connect the extremity to the connector 2 of the LoRa® SPY U.



Case 2: Connection diagram with a Nano SPY U

JRI recommends the use of an IP67 Binder connector (ref 12617) and a universal cable to connect the Nano SPY U to the PT100 Fine probe for Covered Incubators (IVF).

- White wire of the Nano SPY U: Connect to terminal 1 of the IP67 connector.
- Blue wire of the Nano SPY U: Connect to terminal 2 of the IP67 connector.
- Brown wire of the Nano SPY U: Connect to terminal 4 of the IP67 connector.
- White wire of the probe: Connect the extremity to terminal 4 of the IP67 connector.
- 1st Red wire of the probe: Connect the extremity to terminal 1 of the IP67 connector.
- 2nd Red Probe Lead: Connect the extremity to terminal block 2 of the IP67 connector.



Technical drawing

